




WSCD2 rabbit pAb

Catalog No	YP-Ab-08016
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	WSCD2 KIAA0789
Immunogen	Synthesized peptide derived from human WSCD2 AA range: 253-303
Specificity	This antibody detects endogenous levels of WSCD2 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.130% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	WSCD2;WSCD2 KIAA0789;WSC domain-containing protein 2
Observed Band	60kD
Calculated Molecular Weight	60kD
Cell Pathway	Membrane ; Single-pass membrane protein .
Tissue Specificity	Brain,Urinary bladder,
Function	sequence caution:Contaminating sequence. The N-terminus may be contaminated with vector sequence.,similarity:Belongs to the WSCD family.,similarity:Contains 2 WSC domains.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from 293T cells,
primary antibody was diluted at 1:1000, 4° over night