




MUC21 rabbit pAb

Catalog No	YP-Ab-08021
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	MUC21 C6orf205 UNQ697/PRO1342
Immunogen	Synthesized peptide derived from human MUC21 AA range: 34-84
Specificity	This antibody detects endogenous levels of MUC21 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.135% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MUC21;MUC21 C6orf205 UNQ697/PRO1342;Mucin-21 (MUC-21) (Epiglycanin)
Observed Band	60kD
Calculated Molecular Weight	60kD
Cell Pathway	Cell membrane ; Single-pass type I membrane protein ; Extracellular side .
Tissue Specificity	Expressed in lung, large intestine, thymus, and testis. Expressed in normal and malignant bronchial epithelial cells.
Function	miscellaneous:Could be considered as a marker for lung carcinomas.,PTM:O-glycosylated.,tissue specificity:Expressed in lung, large intestine, thymus, and testis. Expressed in normal and malignant bronchial epithelial cells.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from Jarkat cells,
primary antibody was diluted at 1:1000, 4°over night