




OCSTP rabbit pAb

Catalog No	YP-Ab-08022
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	OCSTAMP C20orf123
Immunogen	Synthesized peptide derived from human OCSTP AA range: 488-538
Specificity	This antibody detects endogenous levels of OCSTP at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.136% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	OCSTP;OCSTAMP C20orf123;Osteoclast stimulatory transmembrane protein (OC-STAMP)
Observed Band	60kD
Calculated Molecular Weight	60kD
Cell Pathway	Membrane ; Multi-pass membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from AD293 cells,
primary antibody was diluted at 1:1000, 4° over night