




RN220 rabbit pAb

Catalog No	YP-Ab-08023
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	RNF220 C1orf164
Immunogen	Synthesized peptide derived from human RN220 AA range: 218-268
Specificity	This antibody detects endogenous levels of RN220 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.137% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RN220;RNF220 C1orf164;E3 ubiquitin-protein ligase RNF220 (EC 6.3.2.-) (RING finger protein 220)
Observed Band	60kD
Calculated Molecular Weight	60kD
Cell Pathway	Cytoplasm .
Tissue Specificity	Brain, Eye, Umbilical cord blood,
Function	similarity: Contains 1 RING-type zinc finger.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from U2OS cells,
primary antibody was diluted at 1:1000, 4° over night