




T151B rabbit pAb

Catalog No	YP-Ab-08024
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	TMEM151B C6orf137 TMEM193
Immunogen	Synthesized peptide derived from human T151B AA range: 366-416
Specificity	This antibody detects endogenous levels of T151B at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.138% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	T151B;TMEM151B C6orf137 TMEM193;Transmembrane protein 151B (Transmembrane protein 193)
Observed Band	60kD
Calculated Molecular Weight	60kD
Cell Pathway	Membrane ; Multi-pass membrane protein .
Tissue Specificity	Brain,
Function	similarity:Belongs to the TMEM151 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from THP-1 cells,
primary antibody was diluted at 1:1000, 4° over night