




TTC34 rabbit pAb

Catalog No	YP-Ab-08025
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	TTC34
Immunogen	Synthesized peptide derived from human TTC34 AA range: 352-402
Specificity	This antibody detects endogenous levels of TTC34 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.139% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TTC34;Tetratricopeptide repeat protein 34 (TPR repeat protein 34)
Observed Band	60kD
Calculated Molecular Weight	60kD
Cell Pathway	
Tissue Specificity	
Function	similarity:Contains 8 TPR repeats.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from U2OS cells,
primary antibody was diluted at 1:1000, 4° over night