




CC017 rabbit pAb

Catalog No	YP-Ab-08027
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	C3orf17
Immunogen	Synthesized peptide derived from human CC017 AA range: 396-446
Specificity	This antibody detects endogenous levels of CC017 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.141% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CC017;C3orf17;Uncharacterized protein C3orf17
Observed Band	60kD
Calculated Molecular Weight	60kD
Cell Pathway	Nucleus . Nucleus, nucleolus .
Tissue Specificity	Brain,Cerebellum,Epithelium,Placenta,Testis,
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from KB cells, primary antibody was diluted at 1:1000, 4° over night