




GGT3 rabbit pAb

Catalog No	YP-Ab-08040
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	GGT3P GGT3
Immunogen	Synthesized peptide derived from human GGT3 AA range: 300-350
Specificity	This antibody detects endogenous levels of GGT3 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.154% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	D2S11;dJ18C9.2;Gamma glutamyltransferase 4;Gamma glutamyltransferase 7; Gamma glutamyltransferase like 3;Gamma glutamyltransferase like 5;Gamma glutamyltranspeptidase 7;Gamma glutamyltranspeptidase like 3;GGT 7;GGT4; GGTL3;GGTL5;GGT3;GGT3P GGT3;Putative gamma-glutamyltranspeptidase 3 (GGT 3) (EC 2.3.2.2) (Gamma-glutamyltransferase 3) (Glutathione hydrolase 3) (EC 3.4.19.13) [Cleaved into: Putative gamma-glutamyltranspeptidase 3 heavy chain;Putative gamma-glutamyltranspeptidase 3 light chain]
Observed Band	60kD
Calculated Molecular Weight	60kD
Cell Pathway	Membrane ; Single-pass type II membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from A431 cells,
primary antibody was diluted at 1:1000, 4° over night