




NDNF rabbit pAb

Catalog No	YP-Ab-08041
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	NDNF C4orf31 UNQ2748/PRO6487
Immunogen	Synthesized peptide derived from human NDNF AA range: 75-125
Specificity	This antibody detects endogenous levels of NDNF at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.155% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	NDNF;NDNF C4orf31 UNQ2748/PRO6487;Protein NDNF (Neuron-derived neurotrophic factor)
Observed Band	60kD
Calculated Molecular Weight	60kD
Cell Pathway	Secreted .
Tissue Specificity	Expressed in neurons along the gonadotropin-releasing hormone (GnRH) expressing neurons migratory route.
Function	function:Promotes matrix assembly and cell adhesiveness.,PTM:O-glycosylated; contains heparan sulfate and chondroitin sulfate.,sequence caution:Translation N-terminally extended.,similarity:Contains 2 fibronectin type-III domains., subunit:Binds heparin and chondroitin sulfate.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from DU145 cells, primary antibody was diluted at 1:1000, 4° over night