




ZN358 rabbit pAb

Catalog No	YP-Ab-08046
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	ZNF358
Immunogen	Synthesized peptide derived from human ZN358 AA range: 436-486
Specificity	This antibody detects endogenous levels of ZN358 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.160% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ZN358;ZNF358;Zinc finger protein 358
Observed Band	60kD
Calculated Molecular Weight	60kD
Cell Pathway	Nucleus .
Tissue Specificity	Eye,Placenta,
Function	function:May be involved in transcriptional regulation.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 9 C2H2-type zinc fingers.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from SW480 cells,
primary antibody was diluted at 1:1000, 4° over night