




ZN675 rabbit pAb

Catalog No	YP-Ab-08048
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	ZNF675 TIZ
Immunogen	Synthesized peptide derived from human ZN675 AA range: 66-116
Specificity	This antibody detects endogenous levels of ZN675 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.162% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ZN675;ZNF675 TIZ;Zinc finger protein 675 (TRAF6-binding zinc finger protein) (TRAF6-inhibitory zinc finger protein)
Observed Band	60kD
Calculated Molecular Weight	60kD
Cell Pathway	Nucleus .
Tissue Specificity	Spleen,Thymus,
Function	function:May be involved in transcriptional regulation. May play a role during osteoclast differentiation by modulating TRAF6 signaling activity., similarity:Contains 1 KRAB domain.,similarity:Contains 15 C2H2-type zinc fingers.,subunit:Interacts with TRAF6.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from U2OS cells,
primary antibody was diluted at 1:1000, 4°over night