




EEPD1 rabbit pAb

Catalog No	YP-Ab-08050
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	EEPD1 KIAA1706
Immunogen	Synthesized peptide derived from human EEPD1 AA range: 140-190
Specificity	This antibody detects endogenous levels of EEPD1 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.164% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EEPD1; HSPC107; KIAA1706; EEPD1 KIAA1706; Endonuclease/exonuclease/phosphatase family domain-containing protein 1
Observed Band	60kD
Calculated Molecular Weight	60kD
Cell Pathway	
Tissue Specificity	Blood, Brain, Cord blood, Mammary gland,
Function	similarity: Contains 1 HhH domain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from 3T3 cells, primary antibody was diluted at 1:1000, 4° over night