




NPRL3 rabbit pAb

Catalog No	YP-Ab-08056
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	NPRL3 C16orf35 CGTHBA MARE
Immunogen	Synthesized peptide derived from human NPRL3 AA range: 461-511
Specificity	This antibody detects endogenous levels of NPRL3 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.170% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C16orf35;14 gene protein;-14 gene protein;Alpha-globin regulatory element-containing gene protein;CGTHBA;NPRL3;NPRL3 C16orf35 CGTHBA MARE; Nitrogen permease regulator 3-like protein (-14 gene protein) (Alpha-globin regulatory element-containing gene protein) (Protein CGTHBA)
Observed Band	60kD
Calculated Molecular Weight	60kD
Cell Pathway	Lysosome membrane . Localization to lysosomes is amino acid-independent. .
Tissue Specificity	Widely expressed. Expressed in the frontal lobe cortex as well as in the temporal, parietal, and occipital lobes (PubMed:27173016, PubMed:26505888).
Function	caution:Ser-489 is missing in the human genome assembly but is present in all available mRNAs and ESTs.,similarity:Belongs to the UPF0171 family.,tissue specificity:Widely expressed.,
Background	The function of the encoded protein is not known. [provided by RefSeq, Aug 2011],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from 3T3 cells, primary antibody was diluted at 1:1000, 4° over night