




RGAG4 rabbit pAb

Catalog No	YP-Ab-08057
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	RGAG4 KIAA2001
Immunogen	Synthesized peptide derived from human RGAG4 AA range: 473-523
Specificity	This antibody detects endogenous levels of RGAG4 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.171% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RGAG4;RGAG4 KIAA2001;Retrotransposon gag domain-containing protein 4
Observed Band	60kD
Calculated Molecular Weight	60kD
Cell Pathway	
Tissue Specificity	Amygdala,Brain,
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from 293T cells,
primary antibody was diluted at 1:1000, 4° over night