




TTC22 rabbit pAb

Catalog No	YP-Ab-08058
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	TTC22
Immunogen	Synthesized peptide derived from human TTC22 AA range: 390-440
Specificity	This antibody detects endogenous levels of TTC22 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.173% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TTC22;Tetratricopeptide repeat protein 22 (TPR repeat protein 22)
Observed Band	60kD
Calculated Molecular Weight	60kD
Cell Pathway	
Tissue Specificity	Carcinoma,
Function	similarity:Contains 7 TPR repeats.,
Background	This gene encodes a protein with seven tetratricopeptide (TPR) repeats. Tetratricopeptide repeat containing motifs are found in a variety of proteins and may mediate protein-protein interactions and chaperone activity. Alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, Jun 2011],
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from SH-SY5Y cells, primary antibody was diluted at 1:1000, 4° over night