




CQ066 rabbit pAb

Catalog No	YP-Ab-08063
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	C17orf66
Immunogen	Synthesized peptide derived from human CQ066 AA range: 505-555
Specificity	This antibody detects endogenous levels of CQ066 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.178% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CQ066;C17orf66;Uncharacterized protein C17orf66
Observed Band	65kD
Calculated Molecular Weight	65kD
Cell Pathway	
Tissue Specificity	Brain,Testis,
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HuvEc cells,
primary antibody was diluted at 1:1000, 4° over night