




F227A rabbit pAb

Catalog No	YP-Ab-08065
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	FAM227A
Immunogen	Synthesized peptide derived from human F227A AA range: 429-479
Specificity	This antibody detects endogenous levels of F227A at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.180% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	F227A;FAM227A;Protein FAM227A
Observed Band	65kD
Calculated Molecular Weight	65kD
Cell Pathway	
Tissue Specificity	Testis,
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from 293T cells,
primary antibody was diluted at 1:1000, 4° over night