



FKBP9 rabbit pAb


Catalog No	YP-Ab-08066
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	FKBP9 FKBP60 FKBP63
Immunogen	Synthesized peptide derived from human FKBP9 AA range: 362-412
Specificity	This antibody detects endogenous levels of FKBP9 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.181% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FK506-binding protein 9; FK506 binding protein 9; EC:5.2.1.8; 63 kDa FKBP; 63 kDa FK506-binding protein; FKBP9; FKBP9 FKBP60 FKBP63; Peptidyl-prolyl cis-trans isomerase FKBP9 (PPIase FKBP9) (EC 5.2.1.8) (63 kDa FK506-binding protein) (63 kDa FKBP) (FKBP-63) (FK506-binding protein 9) (FKBP-9) (Rotamase)
Observed Band	65kD
Calculated Molecular Weight	65kD
Cell Pathway	Endoplasmic reticulum .
Tissue Specificity	Brain, Breast, Kidney, Liver, Lung, PNS, Trachea,
Function	catalytic activity: Peptidylproline (omega=180) = peptidylproline (omega=0)., enzyme regulation: Inhibited by FK506., function: PPIases accelerate the folding of proteins during protein synthesis., PTM: Phosphorylated., similarity: Belongs to the FKBP-type PPIase family., similarity: Contains 2 EF-hand domains., similarity: Contains 4 PPIase FKBP-type domains.,
Background	

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from SW480 cells, primary antibody was diluted at 1:1000, 4° over night