




MSTO1 rabbit pAb

Catalog No	YP-Ab-08070
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	MSTO1 LST005 SLTP005
Immunogen	Synthesized peptide derived from human MSTO1 AA range: 349-399
Specificity	This antibody detects endogenous levels of MSTO1 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.185% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SLTP5;MST;MMYAT;LST5;MSTO1;MSTO1 LST005 SLTP005;Protein misato homolog 1
Observed Band	65kD
Calculated Molecular Weight	65kD
Cell Pathway	Mitochondrion outer membrane . Cytoplasm .
Tissue Specificity	Present in all cell lines tested (at protein level). Widely expressed.
Function	function:Involved in the regulation of mitochondrial distribution and morphology., similarity:Belongs to the misato family.,tissue specificity:Present in all cell lines tested (at protein level). Widely expressed.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HpeG2 cells,
primary antibody was diluted at 1:1000, 4° over night