




ZBT44 rabbit pAb

Catalog No	YP-Ab-08074
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	ZBTB44 BTBD15 ZNF851
Immunogen	Synthesized peptide derived from human ZBT44 AA range: 201-251
Specificity	This antibody detects endogenous levels of ZBT44 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.189% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ZBT44; ZBTB44 BTBD15 ZNF851; Zinc finger and BTB domain-containing protein 44 (BTB/POZ domain-containing protein 15) (Zinc finger protein 851)
Observed Band	65kD
Calculated Molecular Weight	65kD
Cell Pathway	Nucleus .
Tissue Specificity	Blood, Brain, Liver, Testis,
Function	function: May be involved in transcriptional regulation., PTM: Phosphorylated upon DNA damage, probably by ATM or ATR., similarity: Contains 1 BTB (POZ) domain., similarity: Contains 4 C2H2-type zinc fingers.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from Hela cells,
primary antibody was diluted at 1:1000, 4°over night