




SOGA3 rabbit pAb

Catalog No	YP-Ab-08080
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	SOGA3 C6orf174
Immunogen	Synthesized peptide derived from human SOGA3 AA range: 805-855
Specificity	This antibody detects endogenous levels of SOGA3 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.195% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SOGA3;SOGA3 C6orf174;Protein SOGA3
Observed Band	105kD
Calculated Molecular Weight	105kD
Cell Pathway	Membrane ; Single-pass membrane protein .
Tissue Specificity	Brain,PCR rescued clones,Testis,
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HEK293 cells,
primary antibody was diluted at 1:1000, 4° over night