




CA101 rabbit pAb

Catalog No	YP-Ab-08102
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	C1orf101
Immunogen	Synthesized peptide derived from human CA101 AA range: 450-500
Specificity	This antibody detects endogenous levels of CA101 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.217% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CA101;C1orf101;Uncharacterized protein C1orf101
Observed Band	105kD
Calculated Molecular Weight	105kD
Cell Pathway	Cell projection, cilium, flagellum membrane ; Single-pass type I membrane protein . Specifically located in the principal piece of sperm tail. .
Tissue Specificity	Testis,
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from MCF-7 cells,
primary antibody was diluted at 1:1000, 4° over night