




CQ104 rabbit pAb

Catalog No	YP-Ab-08105
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	C17orf104
Immunogen	Synthesized peptide derived from human CQ104 AA range: 218-268
Specificity	This antibody detects endogenous levels of CQ104 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.220% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CQ104;C17orf104;Uncharacterized protein C17orf104
Observed Band	105kD
Calculated Molecular Weight	105kD
Cell Pathway	Cytoplasm . Nucleus . at late pachytene a fraction is nuclear. .
Tissue Specificity	Expressed in fetal ovaries (PubMed:26742488). Expressed in testis (PubMed:28380054).
Function	caution:The sequence shown here is derived from an Ensembl automatic analysis pipeline and should be considered as preliminary data.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from K562 cells,
primary antibody was diluted at 1:1000, 4° over night