




ZN507 rabbit pAb

Catalog No	YP-Ab-08113
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	ZN507 KIAA1084
Immunogen	Synthesized peptide derived from human ZN507 AA range: 236-286
Specificity	This antibody detects endogenous levels of ZN507 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.228% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ZN507;ZNF507 KIAA1084;Zinc finger protein 507
Observed Band	105kD
Calculated Molecular Weight	105kD
Cell Pathway	Nucleus .
Tissue Specificity	Brain,Duodenum,Epithelium,Testis,Uterus,
Function	function:May be involved in transcriptional regulation.,sequence caution:Translation N-terminally extended.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 9 C2H2-type zinc fingers.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from VEC cells,
primary antibody was diluted at 1:1000, 4°over night