




ANO4 rabbit pAb

Catalog No	YP-Ab-08117
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	ANO4 TMEM16D
Immunogen	Synthesized peptide derived from human ANO4 AA range: 898-948
Specificity	This antibody detects endogenous levels of ANO4 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.232% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ANO4;anoctamin 4;TMEM16D;Transmembrane protein 16D;ANO4 TMEM16D; Anoctamin-4 (Transmembrane protein 16D)
Observed Band	105kD
Calculated Molecular Weight	105kD
Cell Pathway	Cell membrane ; Multi-pass membrane protein . Shows an intracellular localization. .
Tissue Specificity	Brain,Prostate,
Function	function:May act as a calcium-activated chloride channel.,similarity:Belongs to the anoctamin family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from Hela cells,
primary antibody was diluted at 1:1000, 4° over night