




K121L rabbit pAb

Catalog No	YP-Ab-08144
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	KIAA1211L C2orf55
Immunogen	Synthesized peptide derived from human K121L AA range: 360-410
Specificity	This antibody detects endogenous levels of K121L at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.259% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	K121L;KIAA1211L C2orf55;Uncharacterized protein KIAA1211-like
Observed Band	105kD
Calculated Molecular Weight	105kD
Cell Pathway	
Tissue Specificity	Brain,Epithelium,
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HEK293 cells,
primary antibody was diluted at 1:1000, 4° over night