



CP062 rabbit pAb


Catalog No	YP-Ab-08146
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	C16orf62 101F10.2
Immunogen	Synthesized peptide derived from human CP062 AA range: 304-354
Specificity	This antibody detects endogenous levels of CP062 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.261% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C16orf62;Chromosome 16 open reading frame 62;CP62_HUMAN; DKFZp313M539;DKFZp434B212;Esophageal cancer associated protein; Esophageal cancer-associated protein;FLJ214;Hypothetical protein LOC572; MGC16824;UPF55 protein C16orf62;CP062;C16orf62 101F10.2;UPF0505 protein C16orf62 (Esophageal cancer-associated protein)
Observed Band	105kD
Calculated Molecular Weight	105kD
Cell Pathway	Membrane ; Single-pass membrane protein . Endosome . Endosome location is dependent of the association with the CCC and WASH complexes. .
Tissue Specificity	Epithelium,Heart,Pancreas,Testis,Tongue,Uterus,
Function	similarity:Belongs to the UPF0505 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HpeG2 cells,
primary antibody was diluted at 1:1000, 4° over night