



INT4 rabbit pAb

Catalog No	YP-Ab-08147
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	INTS4 MSTP093
Immunogen	Synthesized peptide derived from human INT4 AA range: 10-60
Specificity	This antibody detects endogenous levels of INT4 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.262% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	INT4;Proto oncogene Int 4 homolog;Proto oncogene Wnt 3;WNT3;INTS4 MSTP093;Integrator complex subunit 4 (Int4)
Observed Band	105kD
Calculated Molecular Weight	105kD
Cell Pathway	Nucleus .
Tissue Specificity	Aorta,Brain,Muscle,Ovary,Uterus,
Function	function:Component of the Integrator complex, a complex involved in the small nuclear RNAs (snRNA) U1 and U2 transcription and in their 3'-box-dependent processing. The Integrator complex is associated with the C-terminal domain (CTD) of RNA polymerase II largest subunit (POLR2A) and is recruited to the U1 and U2 snRNAs genes.,similarity:Belongs to the Integrator subunit 4 family.,similarity:Contains 8 HEAT repeats.,subunit:Belongs to the multiprotein complex Integrator, at least composed of INTS1, INTS2, INTS3, INTS4, INTS5, INTS6, INTS7, INTS8, INTS9/RC74, INTS10, CPSF3L/INTS11 and INTS12.,
Background	INTS4 is a subunit of the Integrator complex, which associates with the C-terminal domain of RNA polymerase II large subunit (POLR2A; MIM 180660) and mediates 3-prime end processing of small nuclear RNAs U1 (RNU1; MIM



180680) and U2 (RNU2; MIM 180690) (Baillat et al., 2005 [PubMed 16239144])
.[supplied by OMIM, Mar 2008],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from U2OS cells,
primary antibody was diluted at 1:1000, 4° over night