




TBD2B rabbit pAb

Catalog No	YP-Ab-08153
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	TBC1D2B KIAA1055
Immunogen	Synthesized peptide derived from human TBD2B AA range: 755-805
Specificity	This antibody detects endogenous levels of TBD2B at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.268% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TBD2B;TBC1D2B KIAA1055;TBC1 domain family member 2B
Observed Band	105kD
Calculated Molecular Weight	105kD
Cell Pathway	Early endosome .
Tissue Specificity	Brain,Colon,Duodenum,Liver,Testis,
Function	function:May act as a GTPase-activating protein.,similarity:Contains 1 PH domain.,similarity:Contains 1 Rab-GAP TBC domain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from KB cells, primary antibody was diluted at 1:1000, 4° over night