




SKOR1 rabbit pAb

Catalog No	YP-Ab-08157
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	SKOR1 CORL1 FUSSEL15 LBXCOR1
Immunogen	Synthesized peptide derived from human SKOR1 AA range: 322-372
Specificity	This antibody detects endogenous levels of SKOR1 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.272% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SKOR1;SKOR1 CORL1 FUSSEL15 LBXCOR1;SKI family transcriptional corepressor 1 (Functional Smad-suppressing element on chromosome 15) (Fussel-15) (LBX1 corepressor 1) (Ladybird homeobox corepressor 1)
Observed Band	105kD
Calculated Molecular Weight	105kD
Cell Pathway	Nucleus .
Tissue Specificity	Present specifically in cerebellar Purkinje cells (at protein level).
Function	caution:The sequence shown here is derived from an Ensembl automatic analysis pipeline and should be considered as preliminary data.,function:Acts as a transcriptional corepressor of LBX1 (By similarity). Inhibits BMP signaling.,similarity:Belongs to the SKI family.,subunit:Interacts with LBX1 (By similarity). Interacts with SMAD1, SMAD2 and SMAD3.,tissue specificity:Present specifically in cerebellar Purkinje cells (at protein level).,
Background	
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from KB cells, primary antibody was diluted at 1:1000, 4° over night