




ZKSC2 rabbit pAb

Catalog No	YP-Ab-08161
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	ZKSCAN2 ZNF694
Immunogen	Synthesized peptide derived from human ZKSC2 AA range: 812-862
Specificity	This antibody detects endogenous levels of ZKSC2 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.276% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ZKSC2;ZKSCAN2 ZNF694;Zinc finger protein with KRAB and SCAN domains 2 (Zinc finger protein 694)
Observed Band	105kD
Calculated Molecular Weight	105kD
Cell Pathway	Nucleus .
Tissue Specificity	Brain,PCR rescued clones,Thymus,Uterine endothelium,
Function	function:May be involved in transcriptional regulation.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 1 KRAB domain.,similarity:Contains 1 SCAN box domain.,similarity:Contains 6 C2H2-type zinc fingers.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

Western blot analysis of lysates from THP-1 cells,
primary antibody was diluted at 1:1000, 4° over night



ZKSC2