




P12L2 rabbit pAb

Catalog No	YP-Ab-08171
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	POM121L2 POM121L
Immunogen	Synthesized peptide derived from human P12L2 AA range: 555-605
Specificity	This antibody detects endogenous levels of P12L2 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.286% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	P12L2;POM121L2 POM121L;POM121-like protein 2
Observed Band	105kD
Calculated Molecular Weight	105kD
Cell Pathway	nuclear pore,nuclear pore central transport channel,nuclear pore nuclear basket,
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from MCF-7 cells,
primary antibody was diluted at 1:1000, 4° over night