




Z3H7A rabbit pAb

Catalog No	YP-Ab-08172
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	ZC3H7A ZC3H7 ZC3HDC7 HSPC055
Immunogen	Synthesized peptide derived from human Z3H7A AA range: 195-245
Specificity	This antibody detects endogenous levels of Z3H7A at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.287% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Z3H7A;ZC3H7A ZC3H7 ZC3HDC7 HSPC055;Zinc finger CCCH domain-containing protein 7A
Observed Band	105kD
Calculated Molecular Weight	105kD
Cell Pathway	Nucleus .
Tissue Specificity	Epithelium,Liver,Placenta,Skeletal muscle,Skin,Umbilical cord blood,
Function	similarity:Contains 1 C2H2-type zinc finger.,similarity:Contains 3 C3H1-type zinc fingers.,similarity:Contains 3 TPR repeats.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from KB cells, primary antibody was diluted at 1:1000, 4° over night