




EFC13 rabbit pAb

Catalog No	YP-Ab-08176
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	EFCAB13 C17orf57
Immunogen	Synthesized peptide derived from human EFC13 AA range: 245-295
Specificity	This antibody detects endogenous levels of EFC13 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.291% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EFC13;EFCAB13 C17orf57;EF-hand calcium-binding domain-containing protein 13
Observed Band	105kD
Calculated Molecular Weight	105kD
Cell Pathway	nucleus,nucleoplasm,cytoplasm,
Tissue Specificity	Human fetal kidney,Testis,
Function	similarity:Contains 6 EF-hand domains.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from Hela cells,
primary antibody was diluted at 1:1000, 4° over night