



PIWL2 rabbit pAb


Catalog No	YP-Ab-08177
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	PIWIL2 HILI
Immunogen	Synthesized peptide derived from human PIWL2 AA range: 754-804
Specificity	This antibody detects endogenous levels of PIWL2 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.292% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Cancer/testis antigen 8;CT8;Hiwi like;MILI;Miwi like;Piwi like 2;Piwi-like protein 2; Piwil2;PIWL2_HUMAN;PIWL2;PIWIL2 HILI;Piwi-like protein 2 (Cancer/testis antigen 80) (CT80)
Observed Band	105kD
Calculated Molecular Weight	105kD
Cell Pathway	Cytoplasm . Present in chromatoid body. Probable component of the meiotic nuage, also named P granule, a germ-cell-specific organelle required to repress transposon activity during meiosis. .
Tissue Specificity	Expressed in adult testis and in most tumors.
Function	catalytic activity:3-beta-hydroxy-Delta(5)-steroid + NAD(+) = 3-oxo-Delta(5)-steroid + NADH.,catalytic activity:Cholest-5-ene-3-beta,7-alpha-diol + NAD(+) = 7-alpha-hydroxycholest-4-en-3-one + NADH.,disease:Defects in HSD3B7 are the cause of congenital bile acid synthesis defect type 1 (CBAS1) [MIM:607765]; also known as neonatal progressive intrahepatic cholestasis. CBAS1 is due to a primary defect in bile synthesis leading to progressive liver disease. Clinical features include neonatal jaundice, severe intrahepatic cholestasis and cirrhosis., function:Plays a central role during spermatogenesis by repressing transposable elements and prevent their mobilization, which is essential for the germline



integrity. Plays an essential role in meiotic differentiation of spermatocytes, germ cell differentiation and in self-renewal of spermatogonial stem cells. Its presence in oocytes suggests tha

Background	PIWL2 belongs to the Argonaute family of proteins, which function in development and maintenance of germline stem cells (Sasaki et al., 2003 [PubMed 12906857]).[supplied by OMIM, Mar 2008],
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from HeLa cells, primary antibody was diluted at 1:1000, 4° over night