




CHADL rabbit pAb

Catalog No	YP-Ab-08181
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	CHADL SLRR4B
Immunogen	Synthesized peptide derived from human CHADL AA range: 493-543
Specificity	This antibody detects endogenous levels of CHADL at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.296% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CHADL;CHADL SLRR4B;Chondroadherin-like protein
Observed Band	85kD
Calculated Molecular Weight	85kD
Cell Pathway	Secreted . Secreted, extracellular space, extracellular matrix .
Tissue Specificity	Brain,
Function	similarity:Belongs to the small leucine-rich proteoglycan (SLRP) family. Class IV subfamily.,similarity:Contains 1 LRRNT domain.,similarity:Contains 21 LRR (leucine-rich) repeats.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from 3T3 cells, primary antibody was diluted at 1:1000, 4° over night