




PHTF1 rabbit pAb

Catalog No	YP-Ab-08183
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	PHTF1 PHTF
Immunogen	Synthesized peptide derived from human PHTF1 AA range: 122-172
Specificity	This antibody detects endogenous levels of PHTF1 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.298% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PHTF1;PHTF1 PHTF;Putative homeodomain transcription factor 1
Observed Band	85kD
Calculated Molecular Weight	85kD
Cell Pathway	Endoplasmic reticulum membrane ; Multi-pass membrane protein . Golgi apparatus, cis-Golgi network membrane ; Multi-pass membrane protein .
Tissue Specificity	Widely expressed with highest levels in testis.
Function	function:May play a role in transcription regulation.,similarity:Contains 1 basic helix-loop-helix (bHLH) domain.,subunit:Interacts with FEM1B.,tissue specificity:Widely expressed with highest levels in testis.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from MCF-7 cells,
primary antibody was diluted at 1:1000, 4° over night