




ANKF1 rabbit pAb

Catalog No	YP-Ab-08185
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	ANKFN1
Immunogen	Synthesized peptide derived from human ANKF1 AA range: 68-118
Specificity	This antibody detects endogenous levels of ANKF1 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.300% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ANKF1;ANKFN1;Ankyrin repeat and fibronectin type-III domain-containing protein 1
Observed Band	85kD
Calculated Molecular Weight	85kD
Cell Pathway	
Tissue Specificity	Brain,
Function	similarity:Contains 1 fibronectin type-III domain.,similarity:Contains 2 ANK repeats.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from U2OS cells,
primary antibody was diluted at 1:1000, 4° over night