




F16B1 rabbit pAb

Catalog No	YP-Ab-08191
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	FAM160B1 KIAA1600
Immunogen	Synthesized peptide derived from human F16B1 AA range: 211-261
Specificity	This antibody detects endogenous levels of F16B1 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.306% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	F16B1;FAM160B1 KIAA1600;Protein FAM160B1
Observed Band	85kD
Calculated Molecular Weight	85kD
Cell Pathway	
Tissue Specificity	[Isoform 1]: Expressed in all tissues tested, highly expressed brain. ; [Isoform 2]: Only detected at high levels in testis.
Function	similarity:Belongs to the UPF0518 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from U2OS cells,
primary antibody was diluted at 1:1000, 4° over night