




TBC16 rabbit pAb

Catalog No	YP-Ab-08208
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	TBC1D16
Immunogen	Synthesized peptide derived from human TBC16 AA range: 529-579
Specificity	This antibody detects endogenous levels of TBC16 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.323% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TBC16;TBC1D16;TBC1 domain family member 16
Observed Band	85kD
Calculated Molecular Weight	85kD
Cell Pathway	early endosome,cytosol,
Tissue Specificity	Brain,Breast,Skin,Urinary bladder,
Function	function:May act as a GTPase-activating protein for Rab family protein(s)., similarity:Contains 1 Rab-GAP TBC domain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from A431 cells,
primary antibody was diluted at 1:1000, 4° over night