




PDZD4 rabbit pAb

Catalog No	YP-Ab-08214
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	PDZD4 KIAA1444 PDZK4 PDZRN4L
Immunogen	Synthesized peptide derived from human PDZD4 AA range: 4-54
Specificity	This antibody detects endogenous levels of PDZD4 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.329% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KIAA1444;PDZ domain-containing protein 4;PDZ domain-containing RING finger protein 4-like protein;PDZK4;PDZRN4L;PDZD4;PDZD4 KIAA1444 PDZK4 PDZRN4L;PDZ domain-containing protein 4 (PDZ domain-containing RING finger protein 4-like protein)
Observed Band	85kD
Calculated Molecular Weight	85kD
Cell Pathway	Cytoplasm, cell cortex . Mainly localized under the plasma membrane.
Tissue Specificity	Brain-specific. Expressed in fetal and adult brain. Up-regulated in synovial carcinomas.
Function	similarity:Contains 1 PDZ (DHR) domain.,subcellular location:Mainly localized under the plasma membrane.,tissue specificity:Brain-specific. Expressed in fetal and adult brain. Up-regulated in synovial carcinomas.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from U2OS cells, primary antibody was diluted at 1:1000, 4° over night