




TTC14 rabbit pAb

Catalog No	YP-Ab-08223
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	TTC14 KIAA1980 UNQ5813/PRO19630
Immunogen	Synthesized peptide derived from human TTC14 AA range: 76-126
Specificity	This antibody detects endogenous levels of TTC14 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.338% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TTC14;TTC14 KIAA1980 UNQ5813/PRO19630;Tetratricopeptide repeat protein 14 (TPR repeat protein 14)
Observed Band	85kD
Calculated Molecular Weight	85kD
Cell Pathway	
Tissue Specificity	Bone marrow,Brain,Human spinal cord,Placenta,Spleen,
Function	similarity:Belongs to the TTC14 family.,similarity:Contains 1 S1 motif domain.,similarity:Contains 4 TPR repeats.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from U2OS cells,
primary antibody was diluted at 1:1000, 4° over night