



LRFN1 rabbit pAb

Catalog No	YP-Ab-08224
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	LRFN1 KIAA1484 SALM2
Immunogen	Synthesized peptide derived from human LRFN1 AA range: 529-579
Specificity	This antibody detects endogenous levels of LRFN1 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.339% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	LRFN1; LRFN1 KIAA1484 SALM2; Leucine-rich repeat and fibronectin type III domain-containing protein 1 (Synaptic adhesion-like molecule 2)
Observed Band	85kD
Calculated Molecular Weight	85kD
Cell Pathway	Membrane ; Single-pass type I membrane protein . Cell junction, synapse . Cell junction, synapse, postsynaptic density membrane . Detected in excitatory, but not inhibitory, synaptic plasma membrane. .
Tissue Specificity	Brain, Ovary, Uterus,
Function	function: Involved in the regulation of excitatory synapses., PTM: Glycosylated., PTM: Phosphorylated upon DNA damage, probably by ATM or ATR., similarity: Belongs to the LRFN family., similarity: Contains 1 fibronectin type-III domain., similarity: Contains 1 Ig-like (immunoglobulin-like) domain., similarity: Contains 7 LRR (leucine-rich) repeats., subunit: Interacts with DLG1, DLGAP1, AMPA receptors and NMDA receptors (By similarity). Interacts with DLG4.,
Background	




matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from 293T cells, primary antibody was diluted at 1:1000, 4° over night