




TMM8A rabbit pAb

Catalog No	YP-Ab-08227
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	TMEM8A TMEM6 TMEM8
Immunogen	Synthesized peptide derived from human TMM8A AA range: 647-697
Specificity	This antibody detects endogenous levels of TMM8A at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.342% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TMM8A;TMEM8A TMEM6 TMEM8;Transmembrane protein 8A (Protein M83) (Transmembrane protein 6) (Transmembrane protein 8)
Observed Band	85kD
Calculated Molecular Weight	85kD
Cell Pathway	Cell membrane ; Multi-pass membrane protein . Lysosome membrane ; Multi-pass membrane protein .
Tissue Specificity	Expressed in pancreas, placenta, spleen, liver, kidney, bone marrow, peripheral blood leukocytes and tonsil.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from MCF-7 cells,
primary antibody was diluted at 1:1000, 4° over night