




DCST2 rabbit pAb

Catalog No	YP-Ab-08232
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	DCST2
Immunogen	Synthesized peptide derived from human DCST2 AA range: 509-559
Specificity	This antibody detects endogenous levels of DCST2 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.347% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DCST2;DC-STAMP domain-containing protein 2
Observed Band	85kD
Calculated Molecular Weight	85kD
Cell Pathway	Membrane ; Multi-pass membrane protein .
Tissue Specificity	Testis,Thymus,
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from Hela cells,
primary antibody was diluted at 1:1000, 4° over night