




F149A rabbit pAb

Catalog No	YP-Ab-08234
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	FAM149A MST119
Immunogen	Synthesized peptide derived from human F149A AA range: 347-397
Specificity	This antibody detects endogenous levels of F149A at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.349% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	F149A;FAM149A MST119;Protein FAM149A
Observed Band	85kD
Calculated Molecular Weight	85kD
Cell Pathway	
Tissue Specificity	Brain,Heart,
Function	sequence caution:Translation N-terminally extended.,similarity:Belongs to the FAM149 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from MCF-7 cells,
primary antibody was diluted at 1:1000, 4° over night