




LONF1 rabbit pAb

Catalog No	YP-Ab-08235
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	LONRF1 RNF191
Immunogen	Synthesized peptide derived from human LONF1 AA range: 690-740
Specificity	This antibody detects endogenous levels of LONF1 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.350% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	LONF1;LONRF1 RNF191;LON peptidase N-terminal domain and RING finger protein 1 (RING finger protein 191)
Observed Band	85kD
Calculated Molecular Weight	85kD
Cell Pathway	cytosol,
Tissue Specificity	Brain,Hepatoma,Ovary,
Function	similarity:Contains 1 Lon domain.,similarity:Contains 1 RING-type zinc finger.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HEK293 cells,
primary antibody was diluted at 1:1000, 4° over night