




CN037 rabbit pAb

Catalog No	YP-Ab-08239
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	C14orf37
Immunogen	Synthesized peptide derived from human CN037 AA range: 371-421
Specificity	This antibody detects endogenous levels of CN037 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.354% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CN037;C14orf37;Uncharacterized protein C14orf37
Observed Band	85kD
Calculated Molecular Weight	85kD
Cell Pathway	Membrane ; Single-pass type I membrane protein .
Tissue Specificity	Kidney,Neuroblastoma,Placenta,
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from 3T3 cells, primary antibody was diluted at 1:1000, 4° over night