




RN183 rabbit pAb

Catalog No	YP-Ab-08249
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	RNF183
Immunogen	Synthesized peptide derived from human RN183 AA range: 4-54
Specificity	This antibody detects endogenous levels of RN183 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RN183;RNF183
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Endoplasmic reticulum membrane ; Single-pass type IV membrane protein . Endoplasmic reticulum . Golgi apparatus, cis-Golgi network membrane . Lysosome membrane .
Tissue Specificity	Kidney and testis.
Function	similarity:Contains 1 RING-type zinc finger.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from SW480 cells,
primary antibody was diluted at 1:1000, 4° over night