




ZN740 rabbit pAb

Catalog No	YP-Ab-08250
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	ZNF740 TB7
Immunogen	Synthesized peptide derived from human ZN740 AA range: 34-84
Specificity	This antibody detects endogenous levels of ZN740 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MGC6176;OriLyt TD-element-binding protein 7;TB7;Zfp74;Zinc finger protein 74; ZN74_HUMAN;ZNF74;ZN740;ZNF740 TB7
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus .
Tissue Specificity	
Function	function:May be involved in transcriptional regulation.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 3 C2H2-type zinc fingers.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from U2OS cells,
primary antibody was diluted at 1:1000, 4° over night